



## Future State: Bridge Transfusion - Bridge Breast Milk - Browse Bottle History

Cerner Workflow ID: 6975 (v. 10.0) Client Workflow ID: 533

Last updated by Reiser Sanders, Feb 06, 2024 12:25pm (UTC -4 hours)

**Workflow Details:**

Workflow Name: Bridge Transfusion - Bridge Breast Milk - Browse Bottle History

Workflow State: Future State

Workstream: Other

Venue: Acute Care

Client Owner:

Cerner Owner:

Standard: Yes

Related Workflow(s):

Tags:

**Workflow Summary:**

Service Line:

Related Solution(s): Bridge Breast Milk Management

Project Name: Niagara Health System:OPT-0297674:NIAG\_CD Niagara HIS RFP

TestBuilder Script(s):

Cerner Workflow ID: 6975 (v. 10.0)

Client Workflow ID: 533

Workflow Notes:

Introduced By: WS 3.2

Validated By: WS 5.2

**Swim Lane:**

Role(s): Nurse (RN)

Department(s):

Security Position(s): Women's Health - Nurse

Nurse - Neonatology

Women's Health - Lactation Consultant

**Start/Stop [3]**

Description: Access Cerner Bridge

**Decision [14]**

Description: Do you want to view one specific bottle's information?

**Work Step [64]**

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar &amp; clicking the magnifying glass

## Future State: Bridge Transfusion - Bridge Breast Milk - Browse Bottle History

Cerner Workflow ID: 6975 (v. 10.0) Client Workflow ID: 533

Last updated by Reiser Sanders, Feb 06, 2024 12:25pm (UTC -4 hours)

**Work Step [25]**

Description: Click on the bottle number to view fortifier information and division information

**Work Step [7]**

Description: Select Back to return to the Browse Bottle History screen

**Work Step [4]**

Description: Scan the patient's wristband or the storage bin label

**Decision [18]**

Description: Do you want to edit a bottle's information?

**Work Step [12]**

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar & clicking the magnifying glass

**Work Step [29]**

Description: Click on the pencil icon to edit the collection date/time, storage state, and/or dispose information.

**Work Step [22]**

Description: Click Confirm to save changes

**Start/Stop [45]**

Description: Select Back to return to the patient's chart

**Work Step [5]**

Description: Select Browse Bottle History

**Decision [24]**

Description: Do you want to reprint a bottle label?

**Work Step [65]**

Description: Select the desired bottle or search for the bottle by typing the bottle number into the search bar & clicking the magnifying glass

**Work Step [36]**

Description: Click the Print hyperlink

## Future State: Bridge Transfusion - Bridge Breast Milk - Browse Bottle History

Cerner Workflow ID: 6975 (v. 10.0) Client Workflow ID: 533

Last updated by Reiser Sanders, Feb 06, 2024 12:25pm (UTC -4 hours)

**Work Step [38]**

Description: Select the desired printer by scanning its barcode or selecting it from the list

**Decision [40]**

Description: Do you want to view multiple bottles by sorting them?

**Work Step [42]**

Description: Click on the column name that is the desired sort parameter

Comments: E.G. Bottle number, state, expiration date, received date, etc.

**Facilities:**

Facility Name: St. Catharines

Status: Approved

Facility Cerner Owner:

Facility Client Owner: Morosin, Julianne

Authorize Date: Feb 06, 2024

Facility Comments: