



NEPHROLOGISTS

Accessing and Navigating Peritoneal Dialysis PowerPlan

1. Navigate to Table of Contents (Blue Menu) and select **+ Add** next to orders.





2. The Add Order dialogue box opens. In the search bar type in Peritoneal Dialysis and select the appropriate regimen to view. Then select **Done**.



Navigating Peritoneal Dialysis PowerPlan

1. The PowerPlan view opens.

The screenshot shows the PowerPlan interface for a CAPD order. It includes sections for Admission/Transfer/Discharge, Vital Signs/Monitoring, Patient Care, Lines/Tubes/Drains, Dialysis Therapy, and Medications. Each item has a checkbox and a description of the task or instruction.

2. Check or uncheck orders based on your needs.
3. Complete missing required information for each order indicated by the following symbol  by right-click next to the order and clicking .
4. **IMPORTANT:** Ensure correct details for the PD order are entered by completing step 3. Below is the 'Details for CAPD' example:





The screenshot shows the 'Details for CAPD' form with various input fields for patient information and dialysis parameters.

Requested Start Date/Time: [Date/Time] EDT	Ideal Weight (kg): [Input]
*Number of Exchanges: [Dropdown]	*Day Fill Volume (mL): [Input]
*Long Dwell Volume (mL): [Input]	*Long Dwell: [Dropdown]
*Heparin (see Orderset if yes): <input type="radio"/> Yes <input type="radio"/> No	Exchange #1 Volume (mL): [Input]
Exchange #1 Solution Type: [Dropdown]	Exchange #1 Instructions: [Input]
Exchange #2 Volume (mL): [Input]	Exchange #2 Solution Type: [Dropdown]
Exchange #2 Instructions: [Input]	Exchange #3 Volume (mL): [Input]
Exchange #3 Solution Type: [Dropdown]	Exchange #3 Instructions: [Input]
Exchange #4 Volume (mL): [Input]	Exchange #4 Solution Type: [Dropdown]



MANAGING PD POWERPLAN HOSPITAL INFORMATION SYSTEM (HIS)

Signing Peritoneal Dialysis PowerPlan

1. Click  to place the order in a planned state or click  to initiate the plan now.
2. Click Orders for Signature.
3. Click .
4. Click  button to view the ordered Peritoneal Dialysis PowerPlan.