


ADD-ON Test Order Monitor


1. Click on Laboratory Console on App Bar.



2. From the pull down menu, select the Add-On template for the location that you are working at (SCS ADD ON, NFS ADD ON or WS ADD ON).



3. Once the template opens any pending ADD ON orders will populate.




Important Information

You MUST enable the auto refresh setting each time you log on or orders will not populate unless you manually click refresh (see instructions below).

Instructions to enable Auto-Refresh


1. Once the template loads
 - a. Click on Auto Refresh: OFF.
 - b. Ensure Enable box is checked off.
 - c. Enter desired interval for auto-refresh and click on.



The screenshot shows the 'Laboratory Console - Collections TAT' window. At the top right, the 'Auto Refresh' status is 'OFF'. Below this is the 'Collection TAT Monitor' section, which includes a 'Collection TAT List' table and a 'Collection TAT Summary' section. A dialog box titled 'Auto Refresh' is open, showing the 'Enable' checkbox checked, the 'Interval (Minutes)' set to '1', and 'OK' and 'Cancel' buttons.

TAT Location	Collection dt/tm	Status	Priority	Order	Accession	Patient	Nurse Collect	List Information
WS 6E	11/03/2024 19:24	Dispatch	AD	MG	102-24-308-00021	SYSTEMTEST, GLTHREE	No	(None)

2. The Add On collection monitor is now set to auto refresh every minute.



The screenshot shows the 'Laboratory Console - Collections TAT' window after the settings have been updated. The 'Auto Refresh' status is now '1 Minutes'. The 'Collection TAT Monitor' section remains the same as in the previous screenshot.

TAT Location	Collection dt/tm	Status	Priority	Order	Accession	Patient	Nurse Collect	List Information
WS 6E	11/03/2024 19:24	Dispatch	AD	MG	102-24-308-00021	SYSTEMTEST, GLTHREE	No	(None)